


Stop'n GO Double Balloon

Target vertebral body reconstruction with the strength of **FOUR** balloons








Advantages of 2 x 2 = 4 independent balloons

- Each balloon can be controlled separately, with individual pressure levels and volumes
- Controlled 'step-by-step' lifting of the endplates
- Optimized vertebral body reconstruction
- Increased safety




X-ray Picture: Inflation of the double balloon catheter; Source: Joline

How the double balloon technique works













Joline GmbH & Co. KG Neue Rottenburger Strasse 50 72379 Hechingen Germany

phone: +49 (0) 7471 9881-0 fax: +49 (0) 7471 9881-222 e-mail: info@joline.de

www.joline.de

Exclusive US-distributor:

Joline Medical LLC 5-C Kidder Road Chelmsford, MA 01824 U.S.A.

phone: 978-245-8800 fax: 978-299-8000 e-mail: info@joline.us www.joline.us Quality first: Joline is EN ISO 13485:2016 certified. Subject to alterations: 53374/21/02_R00

