


Universal tube clamp

**...lighter and safer
for tubings than
metal clamps**



Universal tube clamp
KL 20 with additional function

Accessories



KL 20

Universal tube clamp

Universal tube clamp for medical tubes with a diameter between 2 mm and 10 mm.

This is not a simple 'plastic clamp', but a clamp offering stability that stands up to any comparison with metal clamps.

This quality product guarantees long-lasting and safe use.



Features and advantages

- High-quality, glass fibre reinforced plastic
- Break-resistant
- Secure closing due to interlocking
- Lightweight
- No metal elements, therefore no oxidation
- Suitable for reprocessing



Additional function

- Easy opening of Luer-Lock connectors
- Protects Luer-Lock connectors from damages (caused by metal clamps)



Reprocessing and reuse

- The universal tube clamp can be used several times
- Cleaning can be performed manually or mechanically (cleaning and disinfection unit)
- Reprocessing instructions enclosed

Article no.	DM
KL 20-50	50 units (non-sterile in plastic bag)



Distributed by:

Joline GmbH & Co. KG
Neue Rottenburger Str. 50
72379 Hechingen
Germany

Phone: +49 7471 9881-0
Fax: +49 7471 9881-222
Email: info@joline.de
www.joline.de



Achim Schulz-Lauterbach VMP GmbH
58640 Iserlohn, Germany

Achim Schulz-Lauterbach VMP GmbH is
certified according to DIN EN ISO 9001: 2015



Subject to technical changes
Subject to alterations: 53563/24/10_R00